

Glass DC Nebulizers

SeaSpray

2mL/min uptake	A21-07-USS2
1mL/min uptake	A21-07-USS1
400µL/min uptake	A21-07-USS04

MicroMist

600µL/min uptake	A21-07-UM06
200µL/min uptake	A21-07-UM02
100µL/min uptake	A21-07-UM01
50µL/min uptake	A21-07-UM005



Slurry

4mL/min uptake	A21-07-US6
----------------	------------

Conikal

2mL/min uptake	A21-07-UC2
1mL/min uptake	A21-07-UC1

DC Nebulizers for HF

DuraMist

1mL/min uptake	A21-07-DM1
400µL/min uptake	A21-07-DM04



OpalMist

2mL/min uptake	A21-07-PFA02
600µL/min uptake	A21-07-PFA06
200µL/min uptake	A21-07-PFA02
100µL/min uptake	A21-07-PFA01
50µL/min uptake	A21-07-PFA005



VeeSpray

0.6-3mL/min uptake	A21-07-CV6
--------------------	------------



IsoMist XR


Programmable Temperature Cyclonic Spray Chamber

For ARCOS II SOP

Glass Spray Chamber	KT-1154-XR
Quartz Spray Chamber	KT-1154Q-XR
PFA Spray Chamber	KT-1154P-XR

For Blue SOP

Glass Spray Chamber	KT-1141-XR
Quartz Spray Chamber	KT-1141Q-XR
PFA Spray Chamber	KT-1141P-XR



Spray Chambers

Tracey

50mL Cyclonic	20-809-2665HE
---------------	---------------



Twister

50mL Cyclonic	20-809-2951HE
---------------	---------------



Scott Spray Chamber

20-809-4053

Straight Quartz Adaptor

21-809-3773



Straight Quartz Adapter with Aux Port

21-809-3804

Spray Chambers for HF

Tracey TFE

50mL Cyclonic	20-809-3717
Suitable for use with HF.	
Requires Transfer Tube 21-809-3727.	



Straight HF Resistant Adaptor

21-809-3727

Torches

D-Torch 30-808-3614

Comprises:

Base and Inner Tube	31-808-3605
Quartz Outer Tube	31-808-3611
Injector Ferrule	31-808-2962
Retaining Ring	31-808-3595
Adaptor for D-Torch B19	31-808-3651
Adaptor for D-Torch B13	31-808-3596
Metal Ball Joint Clip	JC-19
Safety Gloves, Medium	SAFEGLOVE-med

Options:

Sheath Gas Adaptor 21-809-4083

**Tapered Injectors for D-Torch**

Quartz Injector 1.0mm	31-808-0009
Quartz Injector 1.2mm	31-808-3643
Quartz Injector 1.5mm	31-808-0007
Quartz Injector 1.8mm	31-808-8146
Quartz Injector 2.0mm	31-808-0008
Quartz Injector 2.4mm	31-808-3541
Quartz Injector 2.5mm	31-808-3645
Alumina Injector 1.0mm	31-808-0835
Alumina Injector 1.5mm	31-808-0013
Alumina Injector 1.8mm	31-808-0015
Alumina Injector 2.0mm	31-808-0027
Alumina Injector 2.4mm	31-808-0385

**Capillary Injectors for D-Torch**

Alumina Injector 1.8mm 31-808-0087

Spectro Torch 31-808-3405

**Injectors for Quartz Torch**

Quartz Injector 1.2mm	31-808-3410
Quartz Injector 1.8mm	31-808-3391
Quartz Injector 2.0mm	31-808-3411
Quartz Injector 3.0mm	31-808-3412
Alumina Injector 1.8mm	31-808-3409
Alumina Injector 2.5mm	31-808-3390



Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1271

**Gaurdian In-Line Filters**

Guardian in-line non-return gas filter 70-803-1362

Guardian in-line sample filter 70-803-1108

**Trident Internal Standard Kit**

Standard Kit 60-808-1185

Standard Kit for HF [60-808-1150](#)



Autosampler Probes

(See page 12 for complete list of probes)

For Cetac ASX-520

Polyimide sheathed Autosampler Probe
0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre Probe
0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre Probe
0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre Probe
1.0mm ID 70-803-0793

