



**GLASS EXPANSION**

Quality By Design

# Guardian Autosampler Probe

- Robust tip design eliminates crushed and damaged tips due to misalignment.
- Drip-resistance prevents cross contamination of samples, especially with oils.
- PEEK shroud ensures critical alignment within the middle of the test tube with built in height adjustment to set the probe depth.
- Built-in particle filtering holds back particulates from blocking the line.
- Completely inert design, all PEEK and PTFE construction.
- Interchangeable UniFit™ sample lines available in various IDs (e.g. 0.3, 0.50, 0.75 & 1.0mm)
- Designed to suit Cetac, Agilent®, PerkinElmer®, Shimadzu®, Aim Lab, and Thermo Scientific™ Autosamplers.



Scan the below QR Code to watch the Guardian Autosampler Probe Performance Comparison video:

