

## Glass DC Nebulizers

**Conikal**  
1mL/min uptake                      A41-1-UC1S

**SeaSpray**  
2mL/min uptake                      A41-1-USS2S

**Slurry**  
4mL/min uptake                      A41-1-US6

**MicroMist**  
600µL/min uptake                      A41-07-UM06S  
200µL/min uptake                      A41-07-UM02  
100µL/min uptake                      A41-07-UM01  
50µL/min uptake                      A41-07-UM005



## IsoMist XR

**Programmable Temperature Cyclonic Spray Chamber**

IsoMist XR Kit                              KT-2749-XR

Requires peristaltic pump

**Comprises:**

Glass Twister Spray Chamber                      21-809-3186  
Assembly (ROL)                              20-809-5176  
Adaptor    21-809-3998  
Adjustable Stand                              21-809-4911  
Metal Ball Joint Clip                              JC-13  
UniFit Drain (PKT 3)                              UFT-16-75



## DC Nebulizers for HF

**OpalMist**  
1mL/min uptake                      A41-1-PFA1S  
600µL/min uptake                      A41-07-PFA06S  
200µL/min uptake                      A41-07-PFA02  
100µL/min uptake                      A41-07-PFA01  
50µL/min uptake                      A41-07-PFA005



**DuraMist**  
1mL/min uptake                      A41-1-DM1S  
400µL/min uptake                      A41-07-DM04S



**VeeSpray**  
0.6 - 3mL/min uptake                      A41-1-CV6



## Spray Chambers

**Tracey**  
50mL Cyclonic                              20-808-8881HE  
With B14 gravity drain



**Twister**  
50mL Cyclonic                              20-809-9119HE  
With B14 gravity drain



**Twister**  
50mL Cyclonic                              20-809-0222HE



## IsoMist XR

**Programmable Temperature Cyclonic Spray Chamber**

IsoMist XR Kit                              KT-1067-XR

**Comprises:**

Glass Twister Spray Chamber  
With B14 drain                              21-809-3270  
Assembly B14 (LOL)                              20-809-5187  
Adaptor    21-809-3998  
Adjustable Stand                              70-803-2798  
Metal Ball Joint Clip                              JC-13  
Drain Device                                      21-808-8214



**Cinnabar**  
20mL Cyclonic                              20-809-0164HE



**HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber**

Hydride Kit                                      KT-1157



**Spray Chambers**

(continued)

**Drain Device** 21-808-8214



**Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch**

Glass Transfer Tube 86mm 21-809-2846



**Spray Chambers for HF**

**Tracey TFE**

50mL Cyclonic 20-809-2508

With B14 gravity drain  
Suitable for use with HF



**HF Drain Device** 21-809-2540



**Mounting bracket for HF Spray Chamber** 21-809-2839



**Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch**

PFA Transfer Tube 82.5mm 31-808-3071



**Torches**

**Semi Demountable Torch (HF)**

30-808-4335

Suitable for ICPS-7500/8100 and ICPE-9000/9800

**Comprises:**

Torch Body 31-808-4341

Injector Adaptor 31-808-4129

Tapered Alumina Injector 1.5mm 31-808-3202

**Torch Options:**

Tapered Quartz Injector 1.0mm 31-808-3423

Tapered Quartz Injector 1.8mm 31-808-4005

Tapered Quartz Injector 2.0mm 31-808-4006

Tapered Quartz Injector 2.5mm 31-808-3302

Tapered Alumina Injector 2.4mm 31-808-3205

Sapphire Injector 2.0mm 31-808-3378



Viton O-rings for FDT/SDT ball joints (PKT 10) 70-V-011

GazFit Connectors (PKT 4) for 4mm OD side arms GAZ-04



**Quartz Torches**


Quartz Torch Standard 1.5mm Injector (Suitable for 7500/8100) 30-808-1179

Quartz Torch for Organics, 1.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0519

Quartz Torch for High TDS, 1.8mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0518

Quartz Torch for High TDS, 2.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0564

Quartz Mini Torch 1.5mm Injector (Suitable for ICPE-9000/9800) 30-807-0534



## RF Coils

Copper	70-900-6100C
Silver	70-900-6100S
Gold	70-900-6100G



**RF Coil Installation Kit** 70-900-6100T



## Common Accessories

### Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1312




### Guardian In-Line Filters

Guardian DC in-line non-return gas filter	70-803-1942
Guardian in-line sample filter	70-803-1108
Eluo Adaptor for Inline Filter	70-803-1160



### Trident Internal Standard Kits

Trident CT Standard Kit	60-703-1179
Standard Kit for HF	60-808-1150



**Magnifier Inspection Tool** 70-803-1923



## Autosampler Probes (See page 20-21 for complete list of probes)

### For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm	
Probe Connecting Line (Red)	70-803-1803
Polyimide sheathed Autosampler Probe	
0.5mm ID	60-808-1186L
PTFE Encapsulated Carbon Fibre Probe	
0.5mm ID	70-803-0784
PTFE Encapsulated Carbon Fibre Probe	
0.75mm ID	70-803-1443
PTFE Encapsulated Carbon Fibre Probe	
1.0mm ID	70-803-0793



### For Shimadzu AS10

PTFE Sheathed Carbon Fibre Probe	
0.75mm ID with UniFit	70-803-1477

